

PIGM antibody

Catalog Number: 113894

General Information

Product name:

PIGM antibody

Immunogen:

Applications:

ELISA, WB, IHC

Specificity:

human, mouse, rat; other species not tested.

Antibody description:

PIGM Rabbit Polyclonal antibody. Positive IHC detected in human testis. Positive WB detected in HeLa cells, PC-3 cells. Observed molecular weight by Western-blot: 49kd

Properties

Preparation:

This antibody was obtained by immunization of PIGM recombinant protein (Accession Number: NM_145167). Purification method: Antigen affinity purified.

Formulation:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage:

Store at -20°C. DO NOT ALIQUOT

Clonality:

Polyclonal

Ig Type:Rabbit IgG


Applications

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200

Validations




HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113894(PIGM antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human testis using Catalog No:113894(PIGM antibody) at dilution of 1:100 (under 10x lens)

PIGM antibody

Catalog Number: 113894



Immunohistochemical of paraffin-embedded human testis using Catalog No:113894(PIGM antibody) at dilution of 1:100 (under 40x lens)