

PIN4 antibody



Catalog Number: 113908

Product name

PIN4 antibody

Recommended Dilution:

WB: 1:200-1:2000

Immunogen

IHC: 1:20-1:200

[Human PIN4 Recombinant protein \(GST tag\)](#)


Validations

Specificity

human, mouse, rat; other species not tested.

Antibody description

PIN4 Rabbit Polyclonal antibody. Positive IHC detected in human liver. Positive WB detected in mouse liver tissue, HeLa cells, mouse brain tissue, mouse kidney tissue. Observed molecular weight by Western-blot: 14kd




Preparation

This antibody was obtained by immunization of PIN4 recombinant protein (Accession Number: BC005234). Purification method: Antigen affinity purified.

mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:113908(PIN4 antibody) at dilution of 1:500

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Immunohistochemical of paraffin-embedded human liver using Catalog No:113908(PIN4 antibody) at dilution of 1:50 (under 10x lens)

Clonality


Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC



Dilutions

Immunohistochemical of paraffin-embedded

PIN4 antibody



Catalog Number: 113908

human liver using Catalog No:113908(PIN4

antibody) at dilution of 1:50 (under 40x lens)