

PLDN antibody



Catalog Number: 113952

Product name

PLDN antibody

Recommended Dilution:

WB: 1:200-1:2000

Immunogen

IP: 1:200-1:2000

[Human PLDN Recombinant protein \(GST tag\)](#)

IHC: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

PLDN Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive IP detected in K-562 cells. Positive WB detected in K-562 cells, HeLa cells, human brain tissue, human heart tissue, mouse thymus tissue. Observed molecular weight by Western-blot: 25kd




Preparation

This antibody was obtained by immunization of PLDN recombinant protein (Accession Number: NM_012388). Purification method: Antigen affinity purified.

K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:113952(PLDN antibody) at dilution of 1:500

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:113952(PLDN antibody) at dilution of 1:50 (under 10x lens)

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions

PLDN antibody



Catalog Number: 113952



IP Result of anti-PLDN (IP:Catalog No:113952, 3ug; Detection:Catalog No:113952 1:500) with K-562 cells lysate 3600ug.