

# PLEKHB1 antibody



Catalog Number: 113960

## Product name

Polyclonal

PLEKHB1 antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

PLEKHB1 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 50 kDa

## Dilutions


Recommended Dilution:

WB: 1:200-1:1000

## Preparation

This antibody was obtained by immunization of PLEKHB1 recombinant protein (Accession Number: NM\_001130035). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:113960(PLEKHB1 Antibody) at dilution of 1:300