

# Plexin B2 antibody

Catalog Number: 113966

## Product name

Plexin B2 antibody

WB: 1:200-1:2000

IP: 1:500-1:5000

## Specificity

Human, Mouse, Rat; other species not tested.

IF: 1:10-1:100

## Antibody description

Plexin B2 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in A549 cells, HEK-293 cells, HeLa cells, HepG2 cells, human brain tissue, MCF7 cells, mouse liver tissue, mouse ovary tissue. Positive IP detected in mouse liver tissue. Observed molecular weight by Western-blot: 80kd, 240kd

## Preparation

This antibody was obtained by immunization of Plexin B2 recombinant protein (Accession Number: BC004542). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, IF, WB, IP


## Dilutions

Recommended Dilution:


## Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:113966(PLXNB2 antibody) at dilution of 1:800



IP Result of anti-PLXNB2 (IP:Catalog No:113966, 4ug; Detection:Catalog No:113966 1:1000) with mouse liver tissue lysate 4000ug.



# Plexin B2 antibody



Catalog Number: 113966

---

Immunofluorescent analysis of HepG2 cells using  
Catalog No:113966(PLXNB2 Antibody) at dilution

of 1:25 and Alexa Fluor 488-conjugated AffiniPure  
Goat Anti-Rabbit IgG(H+L)