

PON1 antibody

Catalog Number: 114058

Product name

PON1 antibody

Immunogen

[Human PON1 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PON1 Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue. Positive WB detected in mouse liver tissue, human blood tissue. Positive IP detected in mouse liver tissue. Observed molecular weight by Western-blot: 42-45 kDa

Preparation

This antibody was obtained by immunization of PON1 recombinant protein (Accession Number: NM_000446). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


IP: 1:200-1:2000

IHC: 1:20-1:200

Validations





mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:114058(PON1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human liver using Catalog No:114058(PON1 antibody) at dilution of 1:50 (under 10x lens)

PON1 antibody

Catalog Number: 114058



Immunohistochemical of paraffin-embedded human liver using Catalog No:114058(PON1 antibody) at dilution of 1:50 (under 40x lens)

IP Result of anti-PON1 (IP:Catalog No:114058, 4ug; Detection:Catalog No:114058 1:500) with mouse liver tissue lysate 4000ug.