

# PRKAR1A antibody



Catalog Number: 114195

---

## Product name

PRKAR1A antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

PRKAR1A Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, mouse cerebellum tissue. Observed molecular weight by Western-blot: 43-45kd

## Preparation

This antibody was obtained by immunization of PRKAR1A recombinant protein (Accession Number: NM\_212471). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114195(PRKAR1A antibody) at dilution of 1:500