

RALB antibody



Catalog Number: 114530

Product name

RALB antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

RALB Rabbit Polyclonal antibody. Positive WB detected in Y79 cells, human brain tissue, human lung tissue, mouse brain tissue, mouse lung tissue. Observed molecular weight by Western-blot: 23 kDa

Preparation

This antibody was obtained by immunization of RALB recombinant protein (Accession Number: NM_002881). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:114530(RALB Antibody) at dilution of 1:600