

RBBP8 antibody



Catalog Number: 114593

Product name

RBBP8 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

RBBP8 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, A431 cells, Jurkat cells, Raji cells. Observed molecular weight by Western-blot: 120kd

Preparation

This antibody was obtained by immunization of RBBP8 recombinant protein (Accession Number: NM_203291). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:114593(RBBP8 antibody) at dilution of 1:300