

RHBDD1 antibody



Catalog Number: 114653

Product name

RHBDD1 antibody

Immunogen

[Human RHBDD1 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RHBDD1 Rabbit Polyclonal antibody. Positive WB detected in human testis tissue, HEK-293 cells, mouse kidney tissue, mouse testis tissue. Observed molecular weight by Western-blot: 36 kDa

Preparation

This antibody was obtained by immunization of RHBDD1 recombinant protein (Accession Number: BC027900). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:114653(RHBDD1 antibody) at dilution of 1:300