

# SCARA5 antibody



Catalog Number: 114985

---

## General Information

### Product name:

SCARA5 antibody

### Immunogen:

### Applications:

ELISA, WB

### Specificity:

Human, Mouse; other species not tested.

### Antibody description:

SCARA5 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 43kd

## Properties

### Preparation:

This antibody was obtained by immunization of Peptide (Accession Number: NM\_173833). Purification method: Antigen affinity purified.

### Formulation:

PBS with 0.02% sodium azide and 50% glycerol

pH 7.3.

### Storage:

Store at -20°C. DO NOT ALIQUOT

### Clonality:

Polyclonal


### Ig Type:Rabbit IgG

## Applications

Recommended Dilution:

WB: 1:200-1:1000

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:114985(SCARA5 antibody) at dilution of 1:300