

SH3BGRL3 antibody



Catalog Number: 115187

Product name

SH3BGRL3 antibody

Immunogen

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SH3BGRL3 Rabbit Polyclonal antibody. Positive IP detected in Mouse brain. Positive WB detected in mouse lung tissue, mouse pancreas tissue. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of SH3BGRL3 recombinant protein (Accession Number: NM_031286). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:

WB: 1:200-1:1000


IP: N/A

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:115187(SH3BGRL3 antibody) at dilution of

IP & WB Mouse brain tissue



IP-PADI2: IP of Mouse brain tissue lysate(5000ug) incubated with 3ug MW of Target band: 75kd

1:100

IP result of anti-PADI2 (Catalog No:113482 for IP and Detection) with mouse brain tissue lysate.