

SHANK1 antibody



Catalog Number: 115263

General Information

Product name:

SHANK1 antibody

Immunogen:

Applications:

ELISA, WB

Specificity:

Human, Mouse, Rat; other species not tested.

Antibody description:

SHANK1 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 250-300 kDa

Properties

Preparation:

This antibody was obtained by immunization of Peptide (Accession Number: NM_016148). Purification method: Antigen affinity purified.

Formulation:

PBS with 0.02% sodium azide and 50% glycerol

pH 7.3.

Storage:

Store at -20°C. DO NOT ALIQUOT

Clonality:

Polyclonal


Ig Type: Rabbit IgG

Applications

Recommended Dilution:

WB: 1:200-1:1000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115263(SHANK1 antibody) at dilution of 1:300