

SLAM/CD150 antibody



Catalog Number: 115294

Product name

SLAM/CD150 antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

Human; other species not tested.

Applications

ELISA, WB

Antibody description

SLAM/CD150 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, human placenta tissue. Observed molecular weight by Western-blot: 70 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Preparation

This antibody was obtained by immunization of SLAM/CD150 recombinant protein (Accession Number: BC012602). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:115294(SLAMF1 Antibody) at dilution of 1:300

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality