

# SNAI2 antibody



Catalog Number: 115442

## Product name

SNAI2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SNAI2 Rabbit Polyclonal antibody. Positive WB detected in mouse pancreas tissue, HepG2 cells. Positive IHC detected in human kidney tissue. Observed molecular weight by Western-blot: 30 kDa

## Preparation

This antibody was obtained by immunization of SNAI2 recombinant protein (Accession Number: NM\_003068). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

## Validations



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with Catalog No:115442(SNAI2 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human kidney using Catalog No:115442(SNAI2 antibody) at dilution of 1:200 (under 10x lens)