

ST6GAL1 antibody



Catalog Number: 115628

Product name

ST6GAL1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ST6GAL1 Rabbit Polyclonal antibody. Positive IP detected in Raji cells. Positive WB detected in U-937 cells, human liver tissue, mouse skeletal muscle tissue, mouse spleen tissue, Raji cells. Observed molecular weight by Western-blot: 50 kDa

Preparation

This antibody was obtained by immunization of ST6GAL1 recombinant protein (Accession Number: XM_017007065). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



U-937 cells were subjected to SDS PAGE followed by western blot with Catalog No:115628(ST6GAL1 Antibody) at dilution of 1:600



IP Result of anti-ST6GAL1 (IP:Catalog No:115628, 4ug; Detection:Catalog No:115628 1:600) with Raji cells lysate 2800ug.