

# ST6GALNAC1 antibody



Catalog Number: 115629

## Product name

Polyclonal

ST6GALNAC1 antibody

## Ig Type

## Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

ST6GALNAC1 Rabbit Polyclonal antibody. Positive WB detected in human stomach tissue. Observed molecular weight by Western-blot: 66-69kd

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of ST6GALNAC1 recombinant protein (Accession Number: BC022462). Purification method: Antigen affinity purified.

## Validations


## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



human stomach tissue were subjected to SDS PAGE followed by western blot with Catalog No:115629(ST6GALNAC1 antibody) at dilution of 1:1200