

# STK10 antibody



Catalog Number: 115710

## Product name

Polyclonal

STK10 antibody

## Ig Type

## Specificity

Rabbit IgG

Human; other species not tested.

## Applications

## Antibody description

ELISA, WB

STK10 Rabbit Polyclonal antibody. Positive WB detected in SKOV-3 cells, A549 cells, COLO 320 cells. Observed molecular weight by Western-blot: 130 kDa

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Preparation

This antibody was obtained by immunization of STK10 recombinant protein (Accession Number: NM\_005990). Purification method: Antigen affinity purified.

## Validations


## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



SKOV-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:115710(STK10 Antibody) at dilution of 1:600