

TCOF1 antibody



Catalog Number: 115914

Product name

TCOF1 antibody

Recommended Dilution:

WB: 1:500-1:5000

Immunogen

[Human TCOF1 Recombinant protein \(GST tag\)](#)

IHC: 1:20-1:200

IF: 1:20-1:200


Specificity

Human; other species not tested.

Validations

Antibody description

TCOF1 Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue. Positive IF detected in MCF-7 cells. Positive WB detected in HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 220 kDa




Preparation

This antibody was obtained by immunization of TCOF1 recombinant protein (Accession Number: NM_001008657). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115914(TCOF1 antibody) at dilution of 1:1000

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:115914(TCOF1 antibody) at dilution of 1:50 (under 10x lens)

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IF

Dilutions

TCOF1 antibody

Catalog Number: 115914



Immunofluorescent analysis of MCF-7 cells, using TCOF1 antibody Catalog No:115914 at 1:50 dilution and FITC-labeled donkey anti-rabbit IgG(green). Blue pseudocolor = DAPI (fluorescent DNA dye).