

TDRD5 antibody

Catalog Number: 115929

Product name

TDRD5 antibody

Immunogen

[Human TDRD5 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

TDRD5 Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Positive WB detected in mouse testis tissue. Positive IHC detected in mouse testis tissue. Observed molecular weight by Western-blot: ~116 kDa

Preparation

This antibody was obtained by immunization of TDRD5 recombinant protein (Accession Number: NM_001199085). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:1000


IP: 1:200-1:1000

IHC: 1:20-1:200

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:115929(TDRD5 Antibody) at dilution of 1:300





+ - Normal Rabbit IgG
- + TDRD5 Antibody

IP Result of anti-TDRD5 (IP:Catalog No:115929, 5ug; Detection:Catalog No:115929 1:300) with mouse testis tissue lysate 2800ug.

TDRD5 antibody

Catalog Number: 115929



Immunohistochemistry of paraffin-embedded mouse testis tissue slide using Catalog No:115929(TDRD5 Antibody) at dilution of 1:50 (under 10x lens)

Immunohistochemistry of paraffin-embedded mouse testis tissue slide using Catalog No:115929(TDRD5 Antibody) at dilution of 1:50 (under 40x lens)