

TFDP2 antibody



Catalog Number: 116007

Product name

TFDP2 antibody

Immunogen

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TFDP2 Rabbit Polyclonal antibody. Positive IP detected in SH-SY5Y cells. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 49kd

Preparation

This antibody was obtained by immunization of TFDP2 recombinant protein (Accession Number: NM_001178138). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116007(TFDP2 antibody) at dilution of 1:400



IP Result of anti-TFDP2 (IP:Catalog No:116007, 4ug; Detection:Catalog No:116007 1:400) with SH-SY5Y cells lysate 1600ug.