

# TFIIE beta antibody



Catalog Number: 116013

## Product name

TFIIE beta antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TFIIE beta Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, A375 cells. Observed molecular weight by Western-blot: 34kd

## Preparation

This antibody was obtained by immunization of TFIIE beta recombinant protein (Accession Number: XM\_024447138). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116013(TFIIE beta antibody) at dilution of 1:800