

# TLR6 antibody

Catalog Number: 116076

---

## Product name

TLR6 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

TLR6 Rabbit Polyclonal antibody. Positive WB detected in RAW 264.7 cells, HT-1080 cells, human brain tissue, Jurkat cells, mouse thymus tissue. Observed molecular weight by Western-blot: 95-110 kDa

## Preparation

This antibody was obtained by immunization of TLR6 recombinant protein (Accession Number: BC111755). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



RAW264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:116076(TLR6 antibody) at dilution of 1:1000