

# TMEM175 antibody



Catalog Number: 116106

## Product name

TMEM175 antibody

## Immunogen

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TMEM175 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, HEK-293 cells, HeLa cells. Observed molecular weight by Western-blot: 56-58kd

## Preparation

This antibody was obtained by immunization of TMEM175 recombinant protein (Accession Number: XM\_017008701). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:116106(TMEM175 antibody) at dilution of 1:500