

Tie2 antibody

Catalog Number: 116123

Product name

Tie2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Tie2 Rabbit Polyclonal antibody. Positive IP detected in mouse lung tissue. Positive WB detected in mouse lung tissue, mouse liver tissue. Observed molecular weight by Western-blot: 136 kDa

Preparation

This antibody was obtained by immunization of Tie2 recombinant protein (Accession Number: NM_000459). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:116123 (TEK antibody) at dilution of 1:500



IP Result of anti-TEK (IP: Catalog No:116123, 4ug; Detection: Catalog No:116123 1:500) with mouse lung tissue lysate 4000ug.