

# TIMM13 antibody



Catalog Number: 116130

## Product name

TIMM13 antibody

## Immunogen

[Human TIMM13 Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

TIMM13 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 11kd

## Preparation

This antibody was obtained by immunization of TIMM13 recombinant protein (Accession Number: NM\_012458). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF

## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:116130(TIMM13 antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells using Catalog No:116130(TIMM13 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)