

TMEM27 antibody



Catalog Number: 116184

Product name

Polyclonal

TMEM27 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

TMEM27 Rabbit Polyclonal antibody. Positive WB detected in mouse pancreas tissue, human spleen tissue, MCF7 cells, mouse kidney tissue. Observed molecular weight by Western-blot: 20 kDa

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of TMEM27 recombinant protein (Accession Number: NM_020665). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

mouse pancreas tissue were subjected to SDS PAGE followed by western blot with Catalog No:116184(TMEM27 antibody) at dilution of 1:500

Clonality