

TOP2B antibody



Catalog Number: 116248

Product name

TOP2B antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

Specificity


Human, Mouse; other species not tested.

IF: 1:10-1:100

Antibody description

TOP2B Rabbit Polyclonal antibody. Positive IHC detected in human cervical cancer tissue. Positive IF detected in HeLa cells, HepG2 cells. Positive IP detected in HepG2 cells. Positive WB detected in HeLa cells, HepG2 cells, K-562 cells. Observed molecular weight by Western-blot: 180kd

Validations




Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001068). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116248(TOP2B antibody) at dilution of 1:300

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF, IP

Immunohistochemical of paraffin-embedded human cervical cancer using Catalog No:116248(TOP2B antibody) at dilution of 1:100 (under 40x lens)

Dilutions

Recommended Dilution:

WB: 1:500-1:5000




TOP2B antibody

Catalog Number: 116248

Immunofluorescent analysis of HeLa cells, using TOP2B antibody Catalog No:116248 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit

IP Result of anti-TOP2B (IP:Catalog No:116248, 4ug; Detection:Catalog No:116248 1:300) with HepG2 cells lysate 3200ug.



IgG (red).