

TOP3A antibody



Catalog Number: 116249

Product name

TOP3A antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TOP3A Rabbit Polyclonal antibody. Positive IP detected in K-562 cells. Positive WB detected in K-562 cells, HL-60 cells. Observed molecular weight by Western-blot: 100-110 kDa

Preparation

This antibody was obtained by immunization of TOP3A recombinant protein (Accession Number: NM_004618). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:116249(TOP3A antibody) at dilution of 1:1000



IP Result of anti-TOP3A (IP:Catalog No:116249, 4ug; Detection:Catalog No:116249 1:500) with K-562 cells lysate 3600ug.