

Trappc9,NIBP antibody



Catalog Number: 116348

Product name

Trappc9,NIBP antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

human, mouse, rat; other species not tested.

Applications

ELISA, WB

Antibody description

Trappc9,NIBP Rabbit Polyclonal antibody. Positive WB detected in Neuro-2a cells, mouse brain tissue, mouse skeletal muscle tissue, SH-SY5Y cells. Observed molecular weight by Western-blot: 128kd,139kd

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations

Preparation

This antibody was obtained by immunization of Trappc9,NIBP recombinant protein (Accession Number: NM_180662). Purification method: Antigen affinity purified.


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



Neuro-2a cells were subjected to SDS PAGE followed by western blot with Catalog No:116348(Trappc9,NIBP antibody) at dilution of 1:300