

# TSNAX antibody



Catalog Number: 116437

## Product name

TSNAX antibody

## Immunogen

[Human TSNAX Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TSNAX Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue. Positive IF detected in MCF-7 cells.

## Preparation

This antibody was obtained by immunization of TSNAX recombinant protein (Accession Number: NM\_005999). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, IF


## Dilutions

Recommended Dilution:


IHC: 1:20-1:200

IF: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:116437(TSNAX antibody) at dilution of 1:50 (under 10x lens)



Immunofluorescent analysis of MCF-7 cells, using TSNAX antibody Catalog No:116437 at 1:50 dilution and FITC-labeled donkey anti-rabbit IgG(green). Blue pseudocolor = DAPI (fluorescent DNA dye).