

UBC13 antibody



Catalog Number: 116507

Product name

UBC13 antibody

Recommended Dilution:

WB: 1:500-1:5000

Immunogen

IP: 1:200-1:2000

[Human UBE2N Recombinant protein \(GST tag\)](#)

IHC: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

UBC13 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive WB detected in mouse thymus tissue, human heart tissue, Jurkat cells, mouse small intestine tissue. Positive IP detected in MCF-7 cells. Observed molecular weight by Western-blot: 17 kDa



Preparation

This antibody was obtained by immunization of UBC13 recombinant protein (Accession Number: NM_003348). Purification method: Antigen affinity purified.

mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:116507(UBE2N antibody) at dilution of 1:1000

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal


Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IP


Dilutions



IP Result of anti-UBE2N (IP:Catalog No:116507, 4ug; Detection:Catalog No:116507 1:600) with MCF-7 cells lysate 1280ug.



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:116507(UBE2N Antibody) at dilution of 1:200 (under 10x lens)



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:116507(UBE2N Antibody) at dilution of 1:200 (under 40x lens)