

# USP33 antibody

Catalog Number: 116606

## Product name

USP33 antibody

## Immunogen

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

USP33 Rabbit Polyclonal antibody. Positive IHC detected in human medulloblastoma tissue, human prostate cancer tissue. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 107 kDa/ 103 kDa/ 93 kDa

## Preparation

This antibody was obtained by immunization of USP33 recombinant protein (Accession Number: NM\_201626). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116606(USP33 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human medulloblastoma using Catalog No:116606(USP33 antibody) at dilution of 1:50 (under 10x lens)