

USP33 antibody



Catalog Number: 116606

Product name

USP33 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

USP33 Rabbit Polyclonal antibody. Positive IHC detected in human medulloblastoma tissue, human prostate cancer tissue. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 107 kDa/ 103 kDa/ 93 kDa

Preparation

This antibody was obtained by immunization of USP33 recombinant protein (Accession Number: NM_201626). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116606(USP33 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human medulloblastoma using Catalog No:116606(USP33 antibody) at dilution of 1:50 (under 10x lens)