

VCAM-1 antibody



Catalog Number: 116725

Product name

ELISA, IHC

VCAM-1 antibody

Dilutions

Specificity

Recommended Dilution:


Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

Antibody description

VCAM-1 Rabbit Polyclonal antibody. Positive IHC detected in human spleen tissue, human kidney tissue, human lung cancer tissue, human placenta tissue, human tonsillitis tissue.

Validations




Preparation

This antibody was obtained by immunization of VCAM-1 recombinant protein (Accession Number: NM_001078). Purification method: Antigen affinity purified.

Immunohistochemistry of paraffin-embedded human spleen tissue slide using Catalog No:116725(VCAM-1 Antibody) at dilution of 1:50 (under 10x lens)

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Immunohistochemistry of paraffin-embedded human spleen tissue slide using Catalog No:116725(VCAM-1 Antibody) at dilution of 1:50 (under 40x lens)

Applications