

VDR antibody



Catalog Number: 116734

Product name

VDR antibody

Immunogen

[Human VDR Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

VDR Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells. Observed molecular weight by Western-blot: 55 kDa

Preparation

This antibody was obtained by immunization of VDR recombinant protein (Accession Number: XM_047429500). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:116734(VDR Antibody) at dilution of 1:1000