

Vegfa antibody

Catalog Number: 116737

Product name

Vegfa antibody

Specificity

Mouse; other species not tested.

Antibody description

Vegfa Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 45 kDa

Preparation

This antibody was obtained by immunization of Vegfa recombinant protein (Accession Number: NM-001287056). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116737(Vegfa Antibody) at dilution of 1:1500