

VEGFR-1/FLT-1 antibody



Catalog Number: 116739

Product name

VEGFR-1/FLT-1 antibody

WB: 1:200-1:2000

IP: 1:200-1:2000

Specificity

Human, Mouse, Rat; other species not tested.

IF: 1:10-1:100

Antibody description

VEGFR-1/FLT-1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, human placenta tissue, mouse lung tissue. Positive IP detected in A549 cells. Positive IF detected in HeLa cells. Positive FC detected in HUVEC cells. Observed molecular weight by Western-blot: 55-60 kDa

Preparation

This antibody was obtained by immunization of VEGFR-1/FLT-1 recombinant protein (Accession Number: NM_001159920). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IP, IF, FC


Dilutions

Recommended Dilution:

Validations




HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116739(VEGFR1, FLT1 antibody) at dilution of



1:800

IP Result of anti-VEGFR1, FLT1 (IP:Catalog No:116739, 5ug; Detection:Catalog No:116739 1:500) with A549 cells lysate 3300ug.



VEGFR-1/FLT-1 antibody



Catalog Number: 116739

Immunofluorescent analysis of HeLa cells using Catalog No:116739(VEGFR1, FLT1 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

1X10⁶ HUVEC cells were stained with 0.2ug VEGFR1, FLT1 antibody (Catalog No:116739, red) and control antibody (blue). Fixed with 4% PFA blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.

