

VEGFR3 antibody



Catalog Number: 116740

Product name

Polyclonal

VEGFR3 antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

VEGFR3 Rabbit Polyclonal antibody. Positive WB detected in mouse spleen tissue, HeLa cells, mouse heart tissue, PC-3 cells. Observed molecular weight by Western-blot: 125 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_002020). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:116740(FLT4 antibody) at dilution of 1:300

Clonality