

VHL antibody



Catalog Number: 116748

Product name

VHL antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

VHL Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IHC detected in human kidney tissue.

Preparation

This antibody was obtained by immunization of VHL recombinant protein (Accession Number: NM_198156). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, IHC


Dilutions

Recommended Dilution:


IHC: 1:20-1:200

IF: 1:10-1:100

Validations



Immunohistochemical of paraffin-embedded human kidney using Catalog No:116748(VHL antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using VHL antibody Catalog No:116748 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).