

**Product name**

VHL antibody

**Specificity**

Human; other species not tested.

**Antibody description**

VHL Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells.

**Preparation**

This antibody was obtained by immunization of VHL recombinant protein (Accession Number: NM\_198156). Purification method: Antigen affinity purified.

**Formulation**

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

**Storage**

Store at -20°C. DO NOT ALIQUOT

**Clonality**

Polyclonal

**Ig Type**

Rabbit IgG


**Applications**

ELISA, IF

**Dilutions**

Recommended Dilution:

IF: 1:10-1:100

**Validations**

Immunofluorescent analysis of HepG2 cells using Catalog No:116749(VHL Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)