

Product name

VHL antibody

Specificity

Human; other species not tested.

Antibody description

VHL Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells.

Preparation

This antibody was obtained by immunization of VHL recombinant protein (Accession Number: NM_198156). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IF

Dilutions

Recommended Dilution:

IF: 1:10-1:100

Validations

Immunofluorescent analysis of HepG2 cells using Catalog No:116749(VHL Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)