

# Vimentin antibody



Catalog Number: 116754

## Product name

Polyclonal

Vimentin antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

Vimentin Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, mouse heart tissue. Observed molecular weight by Western-blot: 57kd

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of Vimentin recombinant protein (Accession Number: XM\_006717500). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116754(VIM antibody) at dilution of 1:1000

## Clonality