

Vinculin antibody



Catalog Number: 116755

Product name

Vinculin antibody

Immunogen

[Human Vinculin Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

Vinculin Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 117 kDa

Preparation

This antibody was obtained by immunization of Vinculin recombinant protein (Accession Number: NM_003373). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:1000-1:10000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116755(Vinculin Antibody) at dilution of 1:2000