

Vinculin antibody



Catalog Number: 116755

Product name

Polyclonal

Vinculin antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse; other species not tested.

Applications

Antibody description

ELISA, WB

Vinculin Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 117 kDa

Dilutions


Recommended Dilution:

WB: 1:1000-1:10000

Preparation

This antibody was obtained by immunization of Vinculin recombinant protein (Accession Number: NM_003373). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116755(Vinculin Antibody) at dilution of 1:2000

Clonality