

Visfatin antibody



Catalog Number: 116758

Product name

Visfatin antibody

Recommended Dilution:

WB: 1:500-1:5000

Specificity

Human, Mouse, Rat, Zebrafish; other species not tested.


IP: 1:500-1:5000

IHC: 1:20-1:200

Antibody description

Visfatin Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor tissue. Positive IP detected in mouse heart tissue. Positive WB detected in mouse heart tissue, HL-60 cells, mouse liver tissue, rat liver tissue, zebrafish tissue. Observed molecular weight by Western-blot: 46-52 kDa

Validations



Preparation

This antibody was obtained by immunization of Visfatin recombinant protein (Accession Number: BC020691). Purification method: Antigen affinity purified.

Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:116758(NAMPT antibody) at dilution of 1:50

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal


Ig Type

Rabbit IgG

Applications


ELISA, IHC, WB, IP

Dilutions



(under 10x lens)

mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:116758(NAMPT antibody) at dilution of 1:1000



IP Result of anti-NAMPT (IP:Catalog No:116758, 3ug; Detection:Catalog No:116758 1:1000) with mouse heart tissue lysate 4000ug.