

VPREB3 antibody



Catalog Number: 116769

Product name

VPREB3 antibody

Specificity

Human; other species not tested.

Antibody description

VPREB3 Rabbit Polyclonal antibody. Positive WB detected in human spleen tissue, human placenta tissue. Observed molecular weight by Western-blot: 46kd

Preparation

This antibody was obtained by immunization of VPREB3 recombinant protein (Accession Number: NM_013378). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:116769(VPREB3 antibody) at dilution of 1:400