

VPS13B antibody



Catalog Number: 116771

Product name

VPS13B antibody

Specificity

Human; other species not tested.

Antibody description

VPS13B Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IHC detected in human liver tissue.

Preparation

This antibody was obtained by immunization of VPS13B recombinant protein (Accession Number: BC132978). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, IHC


Dilutions

Recommended Dilution:


IHC: 1:20-1:200

IF: 1:10-1:100


Validations



Immunofluorescent analysis of HepG2 cells using Catalog No:116771(VPS13B Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog No:116771(VPS13B Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog No:116771(VPS13B Antibody) at dilution of 1:200 (under 40x lens).