

VPS15 antibody

Catalog Number: 116772

Product name

VPS15 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

VPS15 Rabbit Polyclonal antibody. Positive IF detected in Hela cells. Positive WB detected in HepG2 cells. Observed molecular weight by Western-blot: 150kd

Preparation

This antibody was obtained by immunization of VPS15 recombinant protein (Accession Number: NM_014602). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:116772(PIK3R4 antibody) at dilution of 1:500



Immunofluorescent analysis of Hela cells, using PIK3R4 antibody Catalog No:116772 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).