

# ZADH2 antibody

Catalog Number: 116902

## General Information

### Product name:

ZADH2 antibody

### Immunogen:

### Applications:

ELISA, WB, IHC

### Specificity:

human, mouse, rat; other species not tested.

### Antibody description:

ZADH2 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 40kd

## Properties

### Preparation:

This antibody was obtained by immunization of ZADH2 recombinant protein (Accession Number: BC078661). Purification method: Antigen affinity purified.

### Formulation:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

### Storage:

Store at -20°C. DO NOT ALIQUOT

### Clonality:

Polyclonal

**Ig Type:**Rabbit IgG


## Applications

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116902(ZADH2 antibody) at dilution of 1:1000




Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:116902(ZADH2 antibody) at dilution of 1:50 (under 10x lens)

# ZADH2 antibody

Catalog Number: 116902

---



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:116902(ZADH2 antibody) at dilution of 1:50 (under 40x lens)